





Advanced Automation of SPE Methods


Streamlining and Improving Sample
Preparation Procedures for U.S. EPA
Methods

- Does Your Lab Have a Problem?
- Brief Overview of SPE
- Comparison of LLE and SPE
- Evolution of SPE Technology
- Overview of SPE Automation - Benefits
- Optimization of SPE for new Analytical Procedures





Solid Phase Extraction (SPE) is an extraction process whereby an aqueous sample is filtered through a bed of sorbent particles, the analytes of interest are removed from the aqueous sample, and adsorbed **(concentrated)** onto the sorbent. Once concentrated, the analytes are removed from the sorbent by a series of eluting solvent rinses.




- SPE disk particles are typically 5 – 15 micron in size.
- As the water sample passes through the packing, the compounds of interest are retained by several mechanisms: reverse phase, anion exchange, and cation exchange interactions.
- The small particle sizes allow **very fast flow rates**; typically up to 85 - 125 ml/min, without breakthrough.



- Cartridges use 40 micron particles (size). This limits the flow rates to 1 – 10 ml / min. And typically reduces the sample size to 200 ml
- Low flow rates results in long sample prep times.
- Higher flow rates lead to breakthrough.
- Small surface area is not good for dirty samples; plugging will occur.

Particle Size of Cartridges vs. Disk SPE



- Larger particles in Cartridge allows sample breakthrough if flow rates are too fast
- SPE Disks prevent breakthrough even with fast flow rates

Single Sample
Manifold



Multiple Sample Manifold



Automation
System



Liquid-Liquid (LLE) Separatory Funnel Extraction



Adding solvent to
Sep Funnel



Shaking Sep
Funnel for 2 min



Solvent separation
and collection

10 Technician hours
for 20 samples

2 Technicians spend
5 continuous hours
on task

1-Liter water samples
extracted

Comparison of Disk, Cartridge and LLE

Disk

15-25 ml solvent

Filtration

85-125 ml/min

No emulsions

Handles dirty samples

Up to 8 liter samples

20 min / sample

Cartridge

5-15 ml solvent

Filtration

5-10 ml/min

No emulsions

No dirt (clean)

200-500 ml samples

100 min / sample

LLE

200 - 500 ml solvent

Shaking

Emulsions

Limited dirt

2 liter samples

1-2 hours / sample


- Disk SPE Uses **Less Solvent** than a typical LLE Procedure
- Disk SPE is **Faster**; greater sample throughput
- Disk SPE can handle **dirty samples**; river water, waste water, etc
- Disk SPE has **no emulsions**
- Disk SPE can process **large sample volumes** (8 liters)
- Disk SPE Provides **comparable / better recoveries than LLE**
- **New methods** are constantly being developed for SPE

- Increases productivity by improving sample throughput
- Improves accuracy and precision by eliminating variations due to human operators
- Reduces laboratory costs by reducing operator time, glassware, and solvent usage
- Reduces worker exposure to solvents
- Consistent control over critical rinse and air dry steps


- JT Baker Speedisks
- 3M Empore Disk
- Atlantic SPE Disks - *New*
 - Certified for Automation



Pre-Filter on Top
for Dirty Samples



Pace Laboratory



Pace Laboratory

SPE SAMPLE COST ANALYSIS WORKSHEET

This worksheet can be used to help calculate the cost per sample, and allows a cost comparison between a JT Baker Speedisk and a 90 mm CPI C-18 disk. The labor time to set up the extractor, with either disk, would be the same. The only time difference will be in the total time required to filter the sample. This is addressed in the Sample Throughput Worksheet. Clean samples would use the JT Baker Speedisk while dirty samples would use the 90 mm CPI disk. Follow the directions below to fill in the worksheet.

Directions:

- 1) Enter in the cost / ml for each of the reagents to be used.
- 2) If FilterAid400 is to be used, enter 4 under the ML column for the Baker table, and 16 under the ML column for the CPI table.
- 3) If a PreFilter is to be used with the 90 mm disk, enter "1" in the PreFilter line in the 90 mm table.
- 4) Enter the above information in the "green" cells.

Baker C-18 Speedisks	Cost		
REAGENTS	\$	ML	TOTAL
Methylene chloride	0.05	20	1.00
Acetone	0.05	10	0.50
Methanol	0.05	5	0.25
Reagent Water	0.05	5	0.25
			\$ 2.00
DISPOSABLES			
Disk	5.95	1	5.95
FilterAid 400 (4 gms)	0.184	0	0.00
			\$ 5.95
LABOR COST			
Operator time per sample (in minutes)	5		\$ 1.25
TOTAL COST PER SAMPLE			\$ 9.20

CPI 90 mm Disks	Cost		
REAGENTS	\$	ML	TOTAL
Methylene chloride	0.05	80	4.00
Acetone	0.05	40	2.00
Methanol	0.05	20	1.00
Reagent Water	0.05	20	1.00
			\$ 8.00
DISPOSABLES			
Disk	14.79	1	14.79
FilterAid 400 (16 gms)	0.184	0	0.00
PreFilter (GMF-150)	1.75	0	0.00
			\$ 14.79
LABOR COST			
Operator time per sample (in minutes)	5		\$ 1.25
TOTAL COST PER SAMPLE			\$ 24.04

SPE SAMPLE THROUGHPUT and "BREAK-EVEN" WORKSHEET

This spreadsheet is used to calculate the run time per sample, the number of extractors needed for the expected work load, and the estimated break-even time period.

The SPE costs per run, calculated from the Cost Analysis Worksheet will automatically be entered into this Worksheet.

- Directions:**
- 1) Enter in the estimated time to filter a sample. This time will be a function of the sample size; 250 or 1,000 ml.
 - 2) Enter in the necessary Air Dry Filtered Disk time, but do not exceed 10 minutes.
 - 3) Enter in the total number of Extractors, and confirm the desired work load is achieved.
 - 4) See the Purple line for the calculated Breakeven period (in days).
 - 5) Enter in the "Expected" number of samples to be run per day.
 - 6) Enter the above information in the "green" cells

Time Per Sample Calculation			
Prewet Steps (time in sec's)	Dispense	Soak	Dry
Methylene chloride	5	30	30
Acetone	5	30	30
Methanol	5	30	0
Reagent Water	5	30	0
Time (in minutes)	0.33	2.00	1.00
Total time for Prewet Steps			3.33
Filter Sample			10
average estimated time			
Air Dry Filtered Disk			2
max air dry time not to exceed 10 min			
Rinse Steps (time in sec's)			
Acetone	5	60	30
Methylene Chloride	5	60	30
Methylene Chloride	5	60	30
Methylene Chloride	5	60	60
Time (in minutes)	0.33	4	2.5
Total time for Prewet Steps			6.83
TOTAL TIME PER SAMPLE			22.17


Samples Run / Shift / Day Calculation	
1) Enter in number of Horizon extractors	= 2
Calculated time per sample (in minutes) (equal "Total Time per Sample" plus 10 minutes prep)	32
Max number of samples per hour	= 4
Max number of samples per 8 hr. shift	= 30
Max number of samples per 24 hr. day	= 90
Expected number of samples per day	= 30
Price, based on the number of Extractors	= \$25,400
Current Cost to run Sample	= \$ 35.00
	Baker Disk
	50 mm
	CPI Disk
	90 mm
SPE Cost per Run	\$ 9.20
	\$ 24.04
Savings per Sample with SPE	\$ 25.80
	\$ 10.96
Number of Samples to run for Breakeven	984
	2318
Max Number of Samples Run Per Day	90
	90
Number of Days to Breakeven (Max #)	11.0
	25.9
Number of Days to Breakeven (Exp #)	32.8
	77.3

BENEFIT	IMPACT	DESCRIPTION
Reduced solvent use	\$9336	Reduction in solvent use
Waste Disposal	\$3960	Less waste generated
Reduce Labor	\$9408	Reduction of O.T.
Reduce Solvent Contact	unknown	Safety
Reduce Labor	\$7280	Allow Analysts Flexibility
Reduce Glassware	\$15,055	Alleviate replacement costs
Method 525.2	\$116,780	Maintain profit in-house
IDEM contract work	\$150,000	
ERO work	\$10,000	
Reduce ERO contamination	Unknown	Atmospheric Contamination reduced
Future markets	Unknown	Endocrine Disruptors, Personal Products, Pharmaceuticals
→ → → TOTAL \$321,819 → → ¶		

- Since the introduction of disk SPE, there have been significant evolutionary growth
 - **Constant Improvements in the system hardware**
 - **New SPE Packings**
 - **More Analytical Methods for SPE disks**
 - **Improved ability to handle dirty samples**


- Modular Design
- Run up to 8 Extractors
- Run independently
- Run individual methods
- Use 47, 50, and 90-mm SPE disks
- Designed for all SPE Chemistries
- Preprogrammed EPA methods
- Bench top operation
- Handle 4 liter samples – can be run multiple times for greater sample volume





- Gas to deliver solvents
- Vacuum to remove liquids
- All solvents are contained in separate bottles
- Separates waste solvent and waste water
- All sample and solvent pathways are PTFE, PEEK, and Kalrez
- Automatically rinses sample bottle
- Utilizes liquid sensors; keeps disk wet
- Modular design; handle up to 8 extractors.
- Handle up to 4 liter samples; better detection limits.

Envision Platform

Extractor 3 

Admin Help



Current Method: **506-Phth Esters** Purge Method

Abort

Cycle	Step	Procedure
Sample	1	Filtering

Elapsed: **00:07** Remain: **00:00**

Status: Sample processing.







- Developed the Auto Barrel to house the DryDisk membrane
- Designed the Auto Barrel to also work on the 4790 extractors
- Use the Auto Barrel without the membrane
- Allows SPE packing material to be added for cartridge specific applications.
- Packing and frits optimized to provide maximum water sample flow rates, but with highest recoveries



- Experimental Work:
- Used Auto Barrel to hold 5.0 gms of coconut charcoal
- Filter 1 liter of water in 40 to 50 minutes
- Used a glass fiber bottom layer to control the water sample flow rate
- No need to set and monitor the vacuum level – self optimizing
- Recovery values of 50 to 70%

SPE Disks come in similar packing materials as SPE Cartridges for a wide range of applications:

- C8 - Octyl
- C-18 - Octadecyl
- SDB - Styrene Divinyl Benzene
- DVB - Hydrophobic
- DVB - Hydrophilic
- Oil & Grease
- Anion Exchange
- Chelating
- Activated Carbon
- **Mixed Phases – with functional groups**

Alpha Labs,
Westborough, MA

SPE-DEX Automated Extractor System

	Set 1 Ext 1- 5	Set 2 Ext 1- 5	Set 3 Ext 1- 5	Ave All Runs			
Analyte	(ng)	(ng)	(ng)	(ng)	Ave %Rec	SD	%RSD
26DANT	22.99	23.08	23.12	23.06	92	0.832	3.607
24DANT	22.94	23.28	23.62	23.28	93	0.864	3.711
HMX	23.21	23.43	22.00	22.88	92	1.146	5.010
RDX	24.49	24.39	23.42	24.10	96	1.129	4.686
PICRIC	24.57	24.92	25.40	24.96	100	1.034	4.143
1,3,5-TNB	18.64	18.83	18.16	18.54	74	0.920	4.963
1,2-DNB (Surrogate)	23.42	23.54	23.59	23.52	94	0.853	3.626
1,3-DNB	24.46	24.56	23.37	24.13	97	1.289	5.342
NB	22.91	23.04	22.16	22.70	91	1.383	6.092
Tetryl	36.28	35.38	35.13	35.59	142	1.476	4.147

**MN Dept. Agriculture Laboratory
St. Paul, MN**

SPE-DEX Automated Extractor System

Compound	% Rec.	RSD	MDL
EPTC	70	13.65	0.14
Propachlor	96	7.39	0.11
Ethafluralin	85	6.97	0.09
Trifluralin	85	4.70	0.06
Phorate	87	8.48	0.11
Dimethoate	108	4.26	0.07
Simazine	94	2.19	0.03
Propazine	95	1.59	0.02
Terbufos	91	6.34	0.09
Fonofos	88	5.30	0.07
Diazinon	98	3.45	0.05
Triallate	86	4.65	0.06
Methyl Parathion	108	3.50	0.06
Malathion	111	1.34	0.02
Chlorpyrifos	99	2.60	0.04
Pendimethalin	96	2.60	0.04

1 Liter River Water Sample



SPE-DEX Automated
Extractor System

Extracting Chlorinated
Pesticides from 1 Liter
River Water Sample

JT Baker Speedisk

Start Time = 0 min

Note the time
clock

River Water Sample from Taiwan



Automated SPE Disk
Extraction
of
Chlorinated
Pesticides

1 Liter Sample

Start Time = 0 min

← Note the time
clock



Automated SPE Disk
Extraction
of
Chlorinated
Pesticides

500 mL Extracted
Elapsed Time = 6 min

← Note the time
clock



Automated SPE Disk
Extraction
of
Chlorinated
Pesticides

Elapsed Time = 13 min

← Note the time
clock

River Water Sample from Taiwan



SPE-DEX Automated
Disk Extraction
Completed

1Liter River Sample
Extracted in 14 Minutes
Unattended operation

Total Time = 14 min

← Note the time
clock



Automated Extraction and
SPE Disks Extract
Sediment & Particulate
Samples

- Fast
- Precise
- No Emulsions

Total Time = 14 min

← **Note the time
clock**

- Drinking Water (EPA - 500 Series)
- Waste Water (EPA - 600 Series)
- Oil and Grease (EPA - 1664A)
- SW-846 (EPA - 3535A - 8000 Series)
- Endocrine Disruptors
- Hormones
- Food and Beverages
- Crops / Soils

- 8061 Phthalate esters
- 8081 TCLP Organochlorine pesticides
- 8082 PCB's
- 8095 Explosives
- 8141 Organophosphorus pesticides
- 8270 TCLP Semi-volatiles
- 8321 TCLP Phenoxyacid herbicides
- 8330 Nitroaromatics / Nitramines

Prewet:

Solvent	Soak Time (min:sec)	Dry Time (min:sec)
• DCM	1:00	0:30
• Acetone	1:00	0:30
• Water	1:00	0:05
• Water	1:00	0:05

Extract Sample: Typically 5 – 10 min depending on % solids.

- Air dry disk: 0:05

Rinse

Solvent	Soak Time (min:sec)	Dry Time (min:sec)
• Acetone	3:00	0:30
• DCM	1:00	0:30
• DCM	1:00	0:30
• DCM	1:00	1:00

Recovery of Phenols – EPA Method 8270

SPE-DEX Automated Extractor System

Spike Level 5 ppb

Compound	SPE Run			Ave	SD	RSD
	1	2	3			
Phenol*	93	90	87	90	3.12	3.45
2-chloro	90	84	88	87	2.65	3.03
2-methyl	103	92	100	98	5.54	5.64
4-methyl	87	77	80	81	5.43	6.69
2,4-dimethyl	56	49	52	52	3.60	6.87
2,4-dichloro	91	81	83	85	5.54	6.53
2,6-dichloro	94	85	87	89	4.79	5.40
4-chloro-3-methyl	105	104	92	100	6.96	6.94
2,4,6-trichloro	94	84	81	86	6.68	7.74
2,4,5-trichloro	90	81	83	85	4.63	5.45
2,3,4,6-tetrachloro	145	108	89	114	28.5	25.06

Recovery of Poly Aromatic Hydrocarbons (PAH's)

SPE-DEX Automated Extractor System

Compound	RT	Ave Rcvry (n=3)	RSD (n=3)
naphthalene	10.18	87.56	6.02
2-methyl naphthalene	11.87	89.02	7.22
1-methyl naphthalene	12.10	90.48	7.72
acenaphthylene	13.97	87.45	2.80
acenaphthene	14.41	78.41	3.63
fluorene	15.63	89.04	4.03
phenanthrene	17.86	89.50	5.40
anthracene	17.97	85.75	9.56
fluoranthene	20.65	91.65	8.34
pyrene	21.15	88.23	2.59
benz(a)anthracene	24.02	102.46	6.86
chrysene	24.10	109.31	1.93
benzo(b)fluoranthene	26.39	102.54	5.59
benzo(k)fluoranthene	26.44	108.54	2.90
Benzo(a)pyrene	27.02	105.38	10.66
indeno(1,2,3-c,d-) perylene	29.09	96.08	13.74
Benz(a,h)anthracene	29.17	99.92	3.06
benzo(g,h,l)perylene	29.51	92.46	9.28

Data Generated by the US EPA using the SPE-DEX 4790 Automated Extractor System

Compounds Listed in Red are POP's

ng/L	Average	Recovery	% STDEV	MDL
TCMX	35.6	71%	8.5%	8.5
alpha-BHC	42.4	85%	4.3%	4.3
Lindane	43.6	87%	4.0%	4.0
beta-BHC	49.0	98%	4.1%	4.1
Heptachlor	42.1	84%	5.6%	5.6
delta-BHC	41.2	82%	5.0%	5.0
Aldrin	38.6	77%	5.2%	5.2
Hept Epoxide	43.5	87%	4.5%	4.5
gamma-Chlordane	41.3	83%	7.4%	7.4
alpha-Chlordane	42.1	84%	6.7%	6.7
Endosulfan I	41.6	83%	3.8%	3.8
p,p'-DDE	35.2	70%	9.5%	9.5
Dieldrin	40.3	81%	5.4%	5.4
Endrin	40.8	82%	7.1%	7.1
p,p'-DDD	37.9	76%	10.3%	10.3
Endosulfan II	41.9	84%	5.1%	5.1
p,p'-DDT	44.8	90%	11.2%	11.2
Endrin Aldehyde	56.2	112%	10.2%	10.2
Endosulfan Sulfate	43.7	87%	5.5%	5.5
Methoxychlor	51.0	102%	9.2%	9.2
Endrin ketone	41.6	83%	7.5%	7.5
Decachlorobiphenyl	42.6	85%	7.9%	7.9

Northeast Analytical, Inc.

2190 Technology Dr., Schenectady, NY

- Rugged, Fast and Cost effective
- Achieves MDLs (1 PPT (ng/L) for 8-liter samples and 9 PPT (ng/L) for 1-liter samples)
- Reproducible results
- Capable of extracting 1 to 8 liter sample volume
- Minimal solvent usage (30ml/sample)
- No glassware to replace or clean
- Unattended operation
- Limited technician time
- Easy to use



1-Liter Analysis - 1.25 ng/L per Congener Spiked

Congener IUPAC#	Conc ng/L	Rcvry (%)	SD (ng/L)	RSD (%)	MDL (2) (ng/L)
22	1.24	98.92	0.06	5.11	0.20
23	1.25	100.01	0.04	3.46	0.14
24	1.24	98.86	0.04	3.00	0.12
25	1.20	96.21	0.02	1.26	0.05
26	1.32	105.38	0.04	2.93	0.12
27	1.24	98.86	0.04	3.00	0.12
28	1.14	91.46	0.06	4.86	0.17
29	1.29	103.00	0.06	4.55	0.18
30	1.16	92.69	0.04	3.51	0.13
31	1.20	96.40	0.05	4.00	0.15
32	1.21	96.91	0.06	5.21	0.20
33	1.32	105.62	0.04	2.73	0.11
34	1.25	100.20	0.04	3.06	0.12

Data From Northeast Analytical Lab using SPE-DEX 4790 Automated Extractor System

- Endocrine Disruptors
 - Nonyl Phenol and Bisphenol A
 - Estrogens
- Microcystins
- Tri-Butyl Tin
- TPH – Total Petroleum Hydrocarbons
- GRO's and DRO's (Gas and Diesel Range Organics)
- SPE – IR for Off-Shore Oil Platforms

- Endocrine Disruptors – Nonyl Phenol
 - SPE Work done in Japan by Horizon distributor
 - 3M SDB-XD 47mm disk
 - 0.1 ppm spiked with 1 liter sample
 - Prewet 1 = 50% Methanol / 50% Acetone
 - Prewet 2 = Methanol
 - Prewet 3 = Reagent Water
 - Wash 1 = Reagent Water
 - Rinse 1 = 50% Methanol / 50% Acetone
 - HPLC = ODS-3 4.6mm ID, 150mm length
 - Mobile phase = 60% acetonitrile and 40% water, pH 2.6
 - Flow rate = 1.0 ml.min
 - UV at 270 nm

Recovery Values

Run 1 = 89.7%

Run 2 = 86.8%

Run 3 = 90.8%

Run 4 = 92.4%

Run 5 = 90.8%

Run 6 = 90.4%

Run 7 = 91.9%

Run 8 = 92.9%

Ave = 90.7%

RSD = 2.1%

- Endocrine Disruptors – Bisphenol A
 - SPE Work done in Japan by distributor
 - 3M SDB-XD 47mm disk
 - 0.1 ppm spiked with 1 liter sample
 - Prewet 1 = 50% Methanol / 50% Acetone
 - Prewet 2 = Methanol
 - Prewet 3 = Reagent Water
 - Wash 1 = Reagent Water
 - Rinse 1 = 50% Methanol / 50% Acetone
 - HPLC = ODS-3 4.6mm ID, 150mm length
 - Mobile phase = 40% acetonitrile and 60% buffer solution, pH 2.7
 - Flow rate = 1.0 ml.min
 - UV at 270 nm

Recovery Values

Run 1 = 93.8%

Run 2 = 92.1%

Run 3 = 98.2%

Run 4 = 96.1%

Run 5 = 94.8%

Run 6 = 93.4%

Run 7 = 92.6%

Ave = 94.4%

RSD = 2.3%

17- α estradiol, 17- β estradiol, estrone, 17 α -ethinylestradiol

BAKERBOND™ Speedisk Extraction Disk C18 XF

The pH of the 2L water samples was adjusted to 7 using HCl or NaOH (1 M).

CONDITIONING STEP

20 ml Acetone
20 ml Methanol
2 x 10 ml Water

SAMPLE LOADING STEP

WASHING STEP

No washing
Disk dried under vacuum for at least 30 min

ELUTION STEP

(max 1 ml/min)
5 ml Acetone
15 ml Methanol

Compounds	Recovery (%) \pm SD*	Recovery (%) \pm SD**
17- α estradiol	105 \pm 20	107 \pm 22
17- β estradiol	104 \pm 25	103 \pm 27
Estrone	108 \pm 21	107 \pm 18
17 α -ethinylestradiol	102 \pm 21	103 \pm 27

JT Baker Application Note

Effect of SPE Dry Time on Recoveries

Microcystins		Dry Time (min)					
		0	5	10	20	40	80
C18 SD	RR	100	103	85	86	78	39
	YR	101	95	83	82	64	22
	LR	99	103	89	89	70	30
SDB-XD	RR	95		84			71
	YR	100		95			70
	LR	102		99			83

Estrone, 17 β -estradiol, estriol, 16 α -hydroxyestrone, 17 α -ethinylestradiol, mestranol, bestradiol, 17-acetate, β -estradiol 3-sulfate, estone 3-sulfate, roxithromycin, clarithromycin, erythromycin

BAKERBOND *Speedisk* Hydrophobic DVB Extraction Disk

CONDITIONING STEP
15 ml Methanol
15 ml Water

SAMPLE LOADING STEP
(10 ml/min)

WASHING STEP
15 ml water
Air dried for 5 min

ELUTION STEP (max 1 ml/min)
15ml *tert*-butyl methyl ether
15 ml methanol

Analytes	Recovery (%)	RSD (%)*
mestranol	105 \pm 12	11
16 α -hydroxyestrone	74 \pm 15	21
17 β -estradiol	98 \pm 9	9
Estrone	105 \pm 16	15
17 α -ethinylestradiol	83 \pm 6	7
17 β -estradiol acetate	109 \pm 15	14
β -estradiol 3-sulfate	92 \pm 8	9
estone 3-sulfate	95 \pm 12	13
estriol	58 \pm 14	24
clarithromycin	82 \pm 6	7
erythromycin	100 \pm 15	15
roxithromycin	79 \pm 8	11

JT Baker Application Note

Paracetamol, Caffeine, *N,N*-diethyl-3-toluamide, Carbamazepine, Oxazepam, Fluoxetine, Metoprolol, Propranolol, Estrone, 17 β -Estradiol, Clofibric acid, Bezafibrate, Ibuprofen, Diclofenac

BAKERBOND™ SDB-1

CONDITIONING STEP

5 ml hexane
5 ml ethylacetate
10 ml methanol
10 ml water

SAMPLE LOADING STEP

(15 ml/min)

Sample (1 L)

WASHING STEP (5 ml/min)

5 ml water
Dried under nitrogen

ELUTION STEP

30 ml methanol

Compound	Recovery (%)
Paracetamol	60
Caffeine	99
<i>N,N</i> -diethyl-3-toluamide	96
Carbamazepine	100
Oxazepam	65
Fluoxetine	69
Metoprolol	81
Propranolol	68
Estrone	92
17 β -Estradiol	96
Clofibric acid	54
Bezafibrate	55
Ibuprofen	46
Diclofenac	42

JT Baker Application Note

- Adjust 100 ml sample, containing 1-2% methanol to pH 4.5
- Use C18 disk
- Prewet 1 = Ethyl Acetate
- Prewet 2 = Methanol
- Prewet 3 = Reagent Water pH 7.0
- Prewet 4 = Reagent Water pH 4.5
- Filter sample through disk
- Dry SPE Disk in desiccator overnight
- Elute with Ethyl Acetate (acidified with HCL)
- Refrigerate extract overnight – for stability

Otis Evans, Betty Jacobs, and Arnold Cohen - US EPA, Cincinnati, OH
Analyst, January 1991, Vol; 116

<http://www.horizontechinc.com>

Automated Solid Phase Extraction Notes

EPA Method 506	Phthalate and Adipate Esters
EPA Method 508.1	Chlorinated Pesticides, Herbicides and Organohalides
EPA Method 515.2 - Note 1 EPA Method 515.2 - Note 2	Chlorinated Acids
EPA Method 525.2 - Note 1 EPA Method 525.2 - Note 2	Semivolatile Organic Compounds
EPA Method 526	Semivolatile Organic Compounds
EPA Method 548.1 - Note 1 EPA Method 548.1 - Note 2	Endothall
EPA Method 549.1	Diquat and Paraquat
EPA Method 550.1 PAH Note	Polycyclic Aromatic Hydrocarbons (PAHs)
EPA Method 552.1 - Note 1 EPA Method 552.1 - Note 2	Haloacetic Acids and Dalapon
EPA Method 608 - Note 1 EPA Method 608 - Note 2	Chlorinated pesticides and PCB's
EPA Method 625	Bases, Neutrals and Acids
EPA Method 1664A - Note 1 EPA Method 1664A - Note 2	N-Hexane Extractable Material; Oil and Grease
EPA Method 1668A	Chlorinated biophenyl congeners by HRGC/HRMS
EPA Method 8015 - Note 1 EPA Method 8015 - Note 2	Total Petroleum Hydrocarbons (TPH)
EPA Method 8061A	Phthalate Esters
EPA Method 8081A	Automated Solid Phase Extractions of Organochlorine Pesticides from Water Using Certified for Automation Atlantic SPE Disks
EPA Method 8081B	Organochlorine Pesticides

EPA Method 8082A - Note 1	Polychlorinated Biphenyls (PCBs)
EPA Method 8082A - Note 2	NEW Application of Automated Solid Phase Extraction for Determination of Polychlorinated Biphenyls (PCB) in Water
EPA Method 8141B - Note 1 EPA Method 8141B - Note 2	Organophosphorus Compounds
EPA Method 8270D - Note 1 EPA Method 8270D - Note 2 EPA Method 8270D - Dual pH	Semivolatile Organic Compounds
EPA Method 8321B	Solvent-Extractable Nonvolatile Compounds
EPA Method 8330A	Nitroaromatics & Nitramines

Unique Studies

[High Throughput Multi-Residue Screening of Drinking Water using the SPE-DEX and Pegasus GC-TOF MS](#)

[Analysis of Bentazone in River Water](#)

DryVap Notes

[Fast Solvent Exchange DryVap® System Procedure for Pesticides and PCBs](#)

[Impact of Excessive Solvent Washing on Recovery Values](#)

[Neutral and List 1 Pesticides with Horizon Technology DryVap® Concentrator System](#)

[Optimizing the DryVap® System for use with High Boiling Solvent Mixtures](#)

[Recoveries of Organochlorine Pesticides Using MeCl with the DryVap™ System](#)

[Solvent Exchange with DryVap® System for Organochlorine Pesticides with GC-ECD](#)

[Using the Learn Mode to Optimize Recoveries With the DryVap™ System](#)

[Concentrating PCB Congeners 10 Thousand Fold with Horizon DryVap™ Automated Concentrator System](#)

Dry Disk Notes

[Study of DryDisk Background Contamination](#)

[An Improved Solvent Drying Technique for Enhanced Recovery of Organic Compounds](#)

GC Racer Notes

[Fast Separation of TRPH and Diesel Fuel with GC Racer](#)

[5 Minute Analysis of Organochlorine Pesticides with GC Racer](#)

- SPE Disks are a proven technique
- Automation is a cost effective solution to laboratory sample throughput problems
- There are significant advances to automating sample preparation for trace analysis
- Wide Range of Applications developed for SPE disks
- New SPE disk packing materials are being developed for enhanced analyte recoveries



Questions?